



## Hair Straighteners - LABOR PRO

### B254 KERATIN 230

plate straighteners ultrasonic treatment keratin



✉ [bbazar@globelife.com](mailto:bbazar@globelife.com)

☎ +39 0331 1706328

🌐 [www.globelife.com/laborprohair/](http://www.globelife.com/laborprohair/)

### B254 KERATIN 230

plate straighteners ultrasonic treatment keratin

KERATIN plate 230, used at a maximum temperature of 180 ° C, is perfect for a smooth ironed, curls and waves. Used to 230 ° C and combined with ultrasonic function, has all the characteristics to perform the function of sealing the keratin treatments.

- Saucers in Metal Gold.

- Digital display.

- Function ultrasound (helps in the breakdown of molecules to a greater attachment to the Keratin treatment).

- Digital thermostat pre- set at 180 ° C for a canonical use of the plate, but with the possibility of increasing the temperature of up to 230 ° C for treatments to Keratin.

- Silicone cable with swivel 360 degrees.